



Evaluation Support - True Generic - Real Savings - Trusted by NHS since 2018

- + Full support given through all stages of evaluation and trials
- + A True Generic, offering identical performance, handling, specification and sizes
- + Worthwhile savings to be achieved with no compromise

Established & trusted since 2018 within Scotland's largest Health Board

PAHACEL® Fibril Absorbable ORC Haemostat

Item Code	NPC Code	Size
PCF11	ELS935	2.6 CM X 5.1 CM
PCF 12		7.6 CM X 10.2 CM
PCF 13		15.2 CM X 22.9 CM
PCF 14		2.5 CM X 5.1 CM
PCF 15	ELS924	5.1 CM X 10.2 CM
PCF 16	ELS903	10.2 CM X 10.2 CM
PCF 17	ELS923	5 CM X 7.5 CM

PAHACEL® Standard Absorbable ORC Haemostat

Item Code	NPC Code	Size
PCS11	ELS916	1.25CM X 5 CM
PCS12		1.3 CM X 5.1 CM
PCS13		1.5 CM X 5 CM
PCS14		2.5 CM X 2.5 CM
PCS15	ELS917	2.5 CM X 5 CM
PCS16	ELS925	5 CM X 7.5 CM
PCS17	ELS909	5 CM X 35 CM
PCS18		5.1 CM X 7.6 CM
PCS19		5.1 CM X 35.6 CM
PCS20	ELS987	7.5 CM X 10 CM
PCS21	ELS910	10 CM X 20 CM
PCS22		10.2 CM X 20.3 CM
PCS23		15 CM X 23 CM
PCS24	ELS988	12.5 CM X 5 CM
PCS25	ELS912	5 CM X 10 CM

PAHACEL® Pillow Absorbable ORC Haemostat

Item Code	NPC Code	Size
PCE 11		3 CM X 3 CM
PCE 12	ELS85023	5 CM X 5 CM
PCE 13		5 CM X 7.5 CM
PCE 14	ELS85024	2.6 CM X 5.1 CM
PCE 15	ELS85025	5.1 CM X 10.2 CM
PCE 16	ELS85026	10.2 CM X 10.2 CM
PCE 17		2.6 CM X 2.6 CM
PCE 18		7.6 CM X 10.2 CM

PAHACEL® Knit Absorbable ORC Haemostat

Item Code	NPC Code	Size
PCK 11	ELS918	2.6 CM X 2.6 CM
PCK 12	ELS896	7.6 CM X 10.2 CM
PCK 13	ELS943	15.2 CM X 22.9 CM
PCK 14		2.5 CM X 5.1 CM
PCK 15		5.1 CM X 10.2 CM
PCK 16		10.2 CM X 10.2 CM
PCK 17		5 CM X 7.5 CM

ABSORBABLE HAEMOSTAT ORC Products / PROPERTIES

Pahacel® Absorbable Haemostats - Properties

Pahacel® Absorbable Haemostat (Oxidised Regenerated Cellulose) is used adjunctively in surgical procedures to assist in the control of capillary, venous and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective.

Pahacel® Standard, Pahacel® Knit, Pahacel® Fibril and Pahacel® Pillow type haemostats can be cut to size for use in endoscopic procedures. Pahacel® Haemostats are perfectly accepted by the organisms and fully absorbed.

• Materials

Made of Oxidised Regenerated Cellulose. Sterile absorbable knitted fabrics prepared by controlled oxidation of regenerated cellulose. The fabric is white with a pale yellow cast and has a faint, caramel-like aroma; has good tensile strength and can be sutured or cut without fraying.

• Anti-Bacterial Effect

Exhibits some bacteriostatic and bactericidal properties with the effect on growth and multiplication of gram(+) and gram(-) organisms including aerobic, anaerobic bacteria and bacteria shown resistant to antibiotics.

• Bio-Degradability

Absorbed by the body within 7-14 days.

• Haemostatic Effect

Provides rapid haemostatis within 2-3 minutes.

• Sterility

Sterilised by Gamma irradiation technique.

• Bio-Compatibility

ISO-10993 Bio-compatibility results are available.

• Irritation

Does not cause irritation to applied body/tissue.

• Classification

Class III classification according to Medical Device Directives.



STANDARD



KNIT



FIBRIL



PILLOW